



## Scholl Freeze Away Skin Tag remover

### Uses

Clinically proven to remove skin tags using cryotherapy freezing technology in as little as one treatment.

### Dosage

Attach the foam applicator stick by firmly twisting it into the opening on the bottom of the white plastic cap. Keeping canister upright, press the white plastic cap down for 3-5 seconds to saturate the applicator. Remove the foam applicator and the white plastic cap from the shield and wait 15 seconds before applying to the skin tag. Insert the white plastic cap into the skin tag holder. Position the skin tag holder behind the skin tag. Slide the skin tag forward until it reaches the base of the holder. Press down on the plastic cap to make contact with the skin tag, then apply continuously for no longer than 40 seconds.

### Format and Ingredients

Aerosol containing:  
Dimethyl Ether, Propane, Isobutane,

**Category:**  
Skin Care and Antiseptics

**Manufacturer:**  
SCHOLL'S WELLNESS COMPANY  
LTD

**Pack size:** 80ml

**RRP:** £24.99

**Legal Status:** MD