



## Vizulize Advanced Multi-Purpose Contact Lens Solution

### Uses

To cleanse, rinse, disinfect and lubricate all types of contact lenses, including hard, soft and silicone hydrogel lenses.

### Dosage

Wash hands. Place the lens into the palm of the hand, add 3 to 4 drops of solution to each side of the lens surface and gently rub both sides for 20 seconds. Rinse thoroughly with lens solution.

### Format and Ingredients

Solution containing:  
Sodium Hyaluronate, Hypromellose,

**Category:**  
Eye Care

**Manufacturer:**  
EM PHARMA LTD

**Pack size:** 100ml, 360ml

**RRP:** £2.99, £9.99

**Legal Status:** MD, MD