



## Otovent Adult

### Uses

A drug-free treatment for eustachian tube dysfunction.

### Dosage

Inflate the balloon using one nostril while holding the other nostril closed. Use with each nostril up to 3 times a day. Blowing the balloon up creates enough pressure to open the eustachian tube, making it easier for fluid to drain from the middle ear. Use up to 3 times daily until all fluid has drained away.

### Format and Ingredients

Nose piece

**Category:**  
Ear Care

**Manufacturer:**  
KESTREL MEDICAL LTD

**Pack size:** 1 nose piece and  
10 balloons

**RRP:** £14.99

**Legal Status:** MD