



# Endwarts Freeze

## Uses

Uses the cryotherapy method to remove warts and verrucas from hands and feet.

## Dosage

Adults and children over 4 years: once activated, immediately press firmly against the wart for 15 seconds on hands, fingers and arms and for 40 seconds on feet. One treatment is usually sufficient to remove the wart or verruca.

## Format and Ingredients

Freezing agent containing:  
N20 gas,

**Category:**  
Footcare, Skin Care and Antiseptics

**Manufacturer:**  
VIATRIS UK HEALTHCARE LTD

**Pack size:** 7.5g

**RRP:** £19.99

**Legal Status:** MD