



## DENTEK TEMPARIN MAX REPAIR KIT

### Uses

Temporary treatment to replace a lost filling, cap, crown or inlay for up to 48 hours until a dentist is seen.

### Dosage

Follow step-by-step instructions on the pack. A dentist must be visited within 48 hours of use.

### Format and Ingredients

Dental repair kit containing:  
Zinc Oxide, Calcium Sulfate, Paraffinum Liquidum, Potassium Alum, Aluminium Sulfate, Aluminium Phosphate, Menthol, Eugenol, Iron Oxide Yellow,

**Category:**  
Mouth Care

**Manufacturer:**  
PRESTIGE CONSUMER HEALTH

**Pack size:**

**RRP:** £4.99

**Legal Status:** MD